



# VIEW® 3D

3D POST-PROCESSING WORKSTATION

## SYSTEM REQUIREMENTS

	MINIMUM	RECOMMENDED
<b>OS:</b>	Windows 2000, XP, Vista	Windows XP Professional
<b>CPU:</b>	Pentium 1 GHz	Core 2 Quad 2,6 GHz
<b>RAM:</b>	1 GB RAM	2 GB RAM
<b>HDD:</b>	20 GB of empty hard disc space	80 GB of empty hard disc space
<b>Display:</b>	1024x768 pixel	1 or 2 displays with 1280x1024 pixel or more
<b>Graphics:</b>	nVidia GeForce 5600 or ATI X800	nVidia 8800 or ATI Radeon 3870-Serie
<b>LAN:</b>	10 Mbit	100 Mbit
<b>Peripherals:</b>		Scrollmouse CD or DVD writer PostScript printer
<b>Hardware:</b>		Dell hardware

© 2008 IMAGE Information Systems Ltd. - For a Better View in Diagnostics! VIEW® 3D 2.1.0 INT ENG 002R


# THE THIRD DIMENSION



## OUR SOLUTIONS FOR YOUR IMAGING NEEDS

- iQ-VIEW®** The radiology reading station
- iQ-VIEW® 3D** 3D post-processing workstation
- iQ-CAPTURE** Add-on hardware module for capturing images from analog video sources
- iQ-WEBX** DICOM server for storage, teleradiology and image distribution
- iQ-HIS/RIS CLIENT** Simplifying the workflow
- iQ-PRINT** DICOM paper print server
- iQ-ROBOT** Automatic burning and labeling of patient CDs and DVDs
- iQ-ROUTER** Image compression for teleradiology and workflow management
- iQ-WORKLIST** DICOM worklist server optimizing your workflow
- iQ-NUC** Simple nuclear PACS integration
- iQ-RIS** The RIS that simplifies your work
- iQ-DISPLAYS** Medical diagnostic monitors
- iQ-CR ACE** Efficiency in CR
- DICOMReader** Read any DICOM CD into your PACS
- OrthoView™** Add-on module for orthopedic templating and trauma planning





# iQ-VIEW<sup>®</sup> 3D IS A CLINICAL 3D POST-PROCESSING WORKSTATION


... for radiologists. It contains all components of iQ-VIEW<sup>®</sup> with the addition of various CT and MRI image processing features.

Multiplanar Reconstructions (MPR) can be performed at your desktop within seconds and CT and MRI scans can be reformatted and viewed in any desired orientation. Intuitive widgets guide the user to exactly position the viewplane to quickly get the desired viewangle and position in space, therefore no prior training is required.

MPR supports the viewing of trauma, vascular, neurology and oncologic CT and MRI images.

The Maximum Intensity Projection (MIP) maps the densest voxels of CT and MRI scans to your screen. This feature simplifies reading low and high contrast pathologies. MIP is an easy to use feature to display vessels and bones.

Thick Slab Maximum Intensity Projections (thick slab MIP) performs MIPs in a flat layer of CT scans. It is a great feature to assess complex fractures and to easily find even the smallest hair fractures. Never miss a fracture on a trauma CT scan again!



Interactive volume rendering enables the medical user to visualize 3D structures easily, e.g. abdominal vessels, brain stem arteries or organs.

Unwanted structures can be either clipped or cropped individually from the image.

The addition of DENOISE and SOFTEN filters optimize the image impression.

Any visualized images can be easily exported as snapshots or even movies to your PACS, CD/DVD or printer.


Users may perform CT and MRI image post-processing at special post acquisition stations from the device manufacturers, however by using iQ-VIEW<sup>®</sup>3D, you can perform this process faster, easier and keep your CT scanner available by moving your work load from the acquisition consoles to iQ-VIEW<sup>®</sup>3D!

A new and innovative 3D technology makes it possible for the iQ-VIEW<sup>®</sup>3D to run on most standard graphic adapters with a low system resource requirement. iQ-VIEW 3D can even run on a laptop.

iQ-VIEW<sup>®</sup>3D has been proven as the clinical tool of choice by many radiologists worldwide. Ask your dealer for a local reference!

Surface Shaded Display (SSD) visualizes the surface of high contrast objects like bones in 3D e.g. interactively or as rotating sequences.



The epitome of 3D processing is the volume rendering (VRT). Highlight in color any vascular, tumor or bony structure for your patients and referring physicians. Various colors, transparency and light setting options for different contexts make the tool easy to use.



1 iQ-VIEW<sup>®</sup> 3D can be integrated into any PACS. It combines all features of a full 2D reading station and a 3D postprocessing station.

2 You can postprocess any kind of volumetric CT or MRI image data.




3 Reconstructed images can be exported as AVI files, to CD/DVD, to a Windows printer, a DICOM imager or sent back to the PACS.



4 The multiplanar views give the viewing physician a quick overview of the volume. Special projections (e.g. Volume Rendering) may be calculated in the right lower window.

5 A double click on each window magnifies the view immediately (e.g. CT Aorta with Angio preset).

6 Surface Shaded Display is a fast and accurate way to extract and visualize features (e.g. bones/skin).



1 The user interface arranges buttons into functional groups to avoid complexity.

2 Using the clipping tool you can easily adjust your volume of interest by simply dragging handles.

3 Using the cropping tool you may cut out unwanted structures in any plane, e.g. bones or organs overlaying vessels.

4 The thick slab feature makes even the smallest hair fractures visible in multislice CT datasets.

5 In iQ-VIEW<sup>®</sup>3D select the first and last image of a stack and within seconds create a stack of multiplanar or rotated projections. You can store these images to your PACS, imagebox or export as an AVI file.

## iQ-VIEW<sup>®</sup> 3D FEATURES

BASED ON iQ-VIEW<sup>®</sup>/PRO\*

### 3D POST PROCESSING

- MPR – any oblique Multiplanar Reconstruction
- MinIP - Minimum Intensity Projection
- SSD – Surface Shaded Display
- VRT – Online Volume Rendering
- Thick Slab rendering (MIP/VRT/SSD) with variable thickness
- Arbitrary volume cropping
- Image filters

### 3D FEATURES

- Definition of different tissues for volume rendering
- Easy selection of the volume of interest
- 3D zoom/pan and center/window
- 3D-measurements
- 3 different windows for simultaneous image processing
- Improved MPR navigation using widgets
- Support of space navigation devices (3D Connexion)

### GENERAL IMAGE PROCESSING

- Creation of MPR reslice images (double oblique)
- Export of animated sequences (AVI)
- Export of secondary capture images to the local image box, filesystem or PACS
- Thickslab MIP/VRT/SSD export functionality

### GRAPHIC ADAPTERS

- Runs on most standard graphic adapters

### OPERATING SYSTEMS

- Windows 2000, XP, Vista

### LANGUAGES

- English, Finnish, French, German, Italian, Japanese, Polish, Portuguese, Russian, Serbian, Spanish, Turkish

### CERTIFICATION

- CE 0482 and FDA 510(k)



\* See iQ-VIEW<sup>®</sup> for further features.